

flexy[®] Pet Rapid Test

CANINE RABIES ANTIBODY ELISA KIT

This Canine Rabies Antibody ELISA Kit is based on an enzyme-linked immunosorbent assay to detect rabies antibodies in canine serum and plasma, moreover, it can be used in the supervision of rabies vaccine antibody levels in dogs.

SPECIFICATION

Product Code	KE0011
Unit Size	96 tests
Principle	<u>ELISA</u>
Target Species	<u>Canine</u>
Target Samples	<u>Serum, Blood</u>
Incubation Time	60min
Performance	This kit is 100% specific, and sensitivity is >99%
Standard Components:	This kit includes: <ul style="list-style-type: none">- Precoated microtiter plate, 96wells- Negative control- Positive control- Enzyme conjugate- 20xWash buffer- Sample buffer- Stop solution- Kit manual
Storage	2-8°C, do not freeze
Validity	12 months when stored properly

NOTE: this product is only for professional use, not for therapeutic purpose. Confirmation of infection shall be done by veterinarian.



Ring Biotechnology Co., Ltd

E-mail: info@ringbio.com diego@nbgen.com Web: www.ringbio.com

Add: Building 3, Zhongtongtai TechnoPark, No. 11, Kechuang 14th St, Beijing 100176, CHINA

Tel: +86-10-56267496 Technical Support & Service: +86-18600362934